

# Portable Pack Work Instruction

**WI-0117 (3)**

A Service Center Work Instruction  
Security Level: Medium

**Course-0247 (1)**

October 25<sup>th</sup>, 2024

Manager, Service Center-ACP

Effective Date: July 31<sup>st</sup>, 2024

Release Date: See Master Control

## Table of Contents

1. Scope .....	3
2. References .....	3
3. Forms .....	3
4. 300200A Omnistim FX2 Portable.....	<b>Error! Bookmark not defined.</b>
5. 300200A-1 Omnistim FX2 Portable Synchrony .....	3
6. 6000400 FX2 Portable Cycle/Walk .....	4
7. Where Posted.....	6
8. Revision History .....	6

**Security Level: Medium | Confidential Records**

This document is the property of Accelerated Care Plus and must be accounted for. Information herein is considered confidential. Do not reproduce it, reveal it to an unauthorized person or send it outside of Accelerated Care Plus without written authorization from Accelerated Care Plus. This printed document is controlled and only valid through the expiration date indicated on the document.

## 1. Scope

1.1 This work instruction applies to Service Center personnel perform packaging for FX2 portables.

## 2. References

2.1 None.

## 3. Forms

3.1 None.

## 4. 300200A-1 Omnistim FX2 Portable Synchrony

4.1 Begin packing 300200A-1 by ensuring all necessary materials are included:

- 1) FX2 Portable Synchrony Stimulator, (P/N 300200A-1)
- 2) FX2 Portable Synchrony Carry Case.
- 3) 9 Volt Battery, (P/N 74747)
- 4) IPM in envelope.
- 5) Asset Tag Label, \*Always match serial on asset tag to serial listed on the device.
- 6) 14 x 10 x 4 box – Not Pictured.




4.2 Verify 300200A-1 case (2) has blue labeling as seen below.

### Security Level: **Medium** | Confidential Records

This document is the property of Accelerated Care Plus and must be accounted for. Information herein is considered confidential. Do not reproduce it, reveal it to an unauthorized person or send it outside of Accelerated Care Plus without written authorization from Accelerated Care Plus. This printed document is controlled and only valid through the expiration date indicated on the document.



4.3 Package device as seen below:



- 4.3.1 Place items 1,3,4,5 in case (2)
- 4.3.2 Wrap case (2) in bubble wrap and place in box (6), stuffing paper to secure
- 4.3.3 Apply asset tag label (7) to box (6)

## 5. 300400 FX2 Portable Cycle/Walk



**Security Level: Medium | Confidential Records**

This document is the property of Accelerated Care Plus and must be accounted for. Information herein is considered confidential. Do not reproduce it, reveal it to an unauthorized person or send it outside of Accelerated Care Plus without written authorization from Accelerated Care Plus. This printed document is controlled and only valid through the expiration date indicated on the document.

- 5.1 Begin packaging 300400A by ensuring all appropriate materials are included, (1) 300400A device, (2) FX2 Portable case labeled with Green FX2 Cycle/Walk, (3) 9V Battery, (5) Copy of IPM (6) Asset tag.



- 5.2 Place items 1, 3, 4 and 5 into the FX2 Portable case labeled as Omnistim FX2 Cycle/Walk as shown to the right.



- 5.3 Package the case with all included materials into the shipping box 14 x 10 x 4 as shown to the right.



**Security Level: Medium | Confidential Records**

This document is the property of Accelerated Care Plus and must be accounted for. Information herein is considered confidential. Do not reproduce it, reveal it to an unauthorized person or send it outside of Accelerated Care Plus without written authorization from Accelerated Care Plus. This printed document is controlled and only valid through the expiration date indicated on the document.

5.4 Close and seal the box with ACP tape and apply shipping labels as shown to the right.



5.5 Apply asset tag item 6 to the box as shown to the right.



## 6. Where Posted

6.1 This document is posted on Master Control.

## 7. Revision History

7.1 For the complete record of change see Master Control.

### Security Level: **Medium** | Confidential Records

This document is the property of Accelerated Care Plus and must be accounted for. Information herein is considered confidential. Do not reproduce it, reveal it to an unauthorized person or send it outside of Accelerated Care Plus without written authorization from Accelerated Care Plus. This printed document is controlled and only valid through the expiration date indicated on the document.