

MEGAPULSE[®] II DIATHERMY PACKING INSTRUCTIONS



----- The Packaging Materials of the Diathermy box are shown above -----



1. Begin by positioning the diathermy head and arm as seen above



2. Unfold the box and position as seen above.



3. Fold the flaps of the box.



4. Using the provided tape strips, tape the seams, ensuring overlap of each strip, and removing the white tabs after applying each strip.



5. Set the box upright.



6. Roll the MegaPulse II Diathermy unit into the box as shown above.



7. Carefully tip the box with Megapulse II unit until it's lying flat. Fold the flaps and tape seams as shown in step 4.



8. Carefully tip the box back upright. Verify box is securely taped. Apply return service shipping label to the top of the box. Please call the Freight Carrier to schedule equipment to be picked up.

For return instructions, see the yellow or green sheet affixed to the return box.