

OMNISOUND® 3000E PACKING INSTRUCTIONS

**NOTE: the transducers are not interchangeable with other units.
You must return the paired transducers with the unit**



----- Packaging will work for either Omnisound® 3000E -----



1. Remove foam inserts from box.



2. Slide foam inserts over front and rear of the Omnisound® unit as demonstrated above.



3. Place unit and packaging face-up in box.



4. Insert charger, cords and transducers into empty spaces around foam packaging.



5. Tape the top of the box and apply return service shipping label to the top. Please call the Freight Carrier to schedule for the equipment to get picked up.

For return instructions, see the yellow or green sheet affixed to the return box.